

**ON REGISTER**

**ATTIS**  
MEDICAL



**RECOVERY RANGE**

## Excell'ICE Cryotherapy Kit – Knee

### Description :

The Excell'ICE knee kit has been specially designed and engineered to facilitate recovery after surgery, trauma, knee pathology or after sporting activity.

It consists of an anterior compressive cryotherapy package, a posterior cryotherapy package, a pump, a Lycra fabric to protect the skin and an airtight freezer bag.

### Indications:

After knee surgery that requires ice application: prosthesis, osteotomy, ligament surgery, meniscal repair, arthroscopy....

After a knee rehabilitation session. After a sports session.

Size choice: Universal size



## Excell'ICE cryotherapy kit – Ankle

### Description :

The Excell'ICE ankle kit has been specially studied and designed to facilitate recovery after surgery, trauma such as sprain, ankle pathology or after a sporting activity. It consists of a compressive cryotherapy package, a pump, a Lycra fabric to protect the skin and an airtight freezer bag.

### Indications:

**Post-operative ankle surgery. Ankle sprain.**

**After an ankle rehabilitation session. After a sports session.**

Size choice: Universal size



# THE EXCELL'ICE

## Excell'ICE Cryotherapy Kit – Foot

### Description :

The Excell'ICE foot compressive cryotherapy kit has been specially studied and designed to facilitate recovery after surgery or trauma to the forefoot.

It consists of a cold bag that allows you to carry out compressive cryotherapy sessions, a protective lycra fabric, a manual pump and an airtight freezer bag.

### Indications:

Post-operative use of foot surgeries that require support and application of ice: hallux valgus, claw toes, Morton. Metatarsal fracture

Size choice: Universal size



## Excell'ICE Cryotherapy Kit – Elbow

### Description :

Excell'ICE Elbow Compressive Cryotherapy Kit has been specially studied and designed to facilitate recovery after surgery, trauma such as sprain, ankle pathology or after sporting activity.

It consists of a compressive cryotherapy package, a pump, a Lycra fabric to protect the skin and an airtight freezer bag.

### Indications:

Post-operative care for elbow surgery, sports trauma (sprain, rupture, tendinopathy, etc.), after rehabilitation session. Analgesic, anti-edematous and anti-hemorrhagic effects.

Size choice: Universal size



## Should'ICE cryotherapy support vest

### Description :

The Should'ICE shoulder support vest was specially studied and designed to facilitate recovery after shoulder surgery. It consists of a shoulder support vest, a cold bag that allows for cryotherapy sessions and an airtight freezer bag. Easy to use, the vest can be put on and used without external assistance.

### Indications:

Surgical sets for shoulder surgery (prosthesis, rotator cuff, pillar, etc.). Shoulder trauma.

### Bilateral model

### Size choice:

chest circumference < 100 cm: S/M  
chest circumference ≥ 100 cm: L/XL



## KHip'Cool Abdominal/Hip Support Belt + Compressive Cryotherapy Pack

### Description :

The KHip'Cool abdominal belt was specially studied and designed to facilitate recovery after abdominal surgery. It consists of a belt that supports the abdomen, a cold bag that allows cryotherapy sessions, a manual pump and an airtight freezer bag. Versatile, the KHip Cool belt can also be used post-operatively after hip surgery thanks to a specific strap.

### Indications:

Postpartum temporary abdominal wall deficiencies, bruises and fractures of the ribs or sternum... Post-operative use for abdominal surgery or hip surgery.

Size choice: Universal size



**ATTIS**  
M E D I C A L

