


EM REGISTRO

ATTIS
MEDICAL

MEDIZIP



The modern way of
closing wounds:
MEDIZIP!

Close Wounds Without Needing Clamps, Stitches or
anesthesia

MEDIZIP



Introducing the MEDIZIP

The MEDIZIP is an atraumatic wound closure method which was conceived by Prof. Dr. Kaeßmann. Suture channels and staple perforations can be dispensed with when the MEDIZIP is used in surgical procedures.

On each side of the wound a medical zipper will be applied. A plaster, specially produced for this function, is the connection between the skin and the zipper. By closing the zipper, the sides of the wound are merged together to ensure a perfect healing process.



Advantages of the MEDIZIP



Reduction in the length of operations and time spent under anaesthesia. No manual assistance or special instruments required.



Significantly reduced risk of infection due to the elimination of retrograde infection routes e.g. suture channels.




No mechanical irritation in the wound opening.



Painless removal of the MEDIZIP from the skin wound after 7 – 11 days, without the necessity for special instruments.



Short learning phase for new users - implementation possible by surgical support staff.



The MEDIZIP can also be used for the closure of curved wounds



As this is a disposable product, the time required for the cleaning and re-sterilisation of the customary instruments is dispensed with.