

ON REGISTER

**ATTIS**  
MEDICAL



# D<sup>2</sup>HAND

DESIGN AND MANUFACTURE FOR  
**ORTHOPEDIC SURGERY**



**DELPHOS**  
Implants

Embracing Life™

# SYSTEM OVERVIEW


1



2



3



**ONE TRAY**  
**FOUR SYSTEMS**  
**ALL LOCKING**



# 1,2 mm

## REGULAR HAND PLATE SYSTEM



locking holes

0,6 mm & 0,8 mm  
thickness plates

t I  0,6 mm

t I  0,8 mm



### Self-tapping 1,2 mm screw

Parafuso autorroscante 1,2 mm  
Vis autotaraudante 1,2 mm  
Tornillo autorroscante 1,2 mm



### Self-tapping locking 1,2 mm screw


Parafuso de bloqueio autorroscante 1,2 mm  
Vis de blocage autotaraudante 1,2 mm  
Tornillo de bloqueo autorroscante 1,2 mm



### Self-tapping LAG 1,2 mm screw


Parafuso LAG autorroscante 1,2 mm  
Vis LAG autotaraudante 1,2 mm  
Tornillo LAG autorroscante 1,2 mm

atraumatic tip




1,2 mm  
Ø thread

- DTH-1204-265** - 4 mm lenght
- DTH-1205-265** - 5 mm lenght
- DTH-1206-265** - 6 mm lenght
- DTH-1207-265** - 7 mm lenght
- DTH-1208-265** - 8 mm lenght
- DTH-1209-265** - 9 mm lenght
- DTH-1210-265** - 10 mm lenght
- DTH-1211-265** - 11 mm lenght
- DTH-1212-265** - 12 mm lenght
- DTH-1213-265** - 13 mm lenght
- DTH-1214-265** - 14 mm lenght
- DTH-1216-265** - 16 mm lenght
- DTH-1218-265** - 18 mm lenght



1,2 mm  
Ø thread

- DTH-1204-265K** - 4 mm lenght
- DTH-1205-265K** - 5 mm lenght
- DTH-1206-265K** - 6 mm lenght
- DTH-1207-265K** - 7 mm lenght
- DTH-1208-265K** - 8 mm lenght
- DTH-1209-265K** - 9 mm lenght
- DTH-1210-265K** - 10 mm lenght
- DTH-1211-265K** - 11 mm lenght
- DTH-1212-265K** - 12 mm lenght
- DTH-1213-265K** - 13 mm lenght
- DTH-1214-265K** - 14 mm lenght
- DTH-1216-265K** - 16 mm lenght
- DTH-1218-265K** - 18 mm lenght




1,2 mm  
Ø thread

- DTH-1206-365** - 6 mm lenght
- DTH-1207-365** - 7 mm lenght
- DTH-1208-365** - 8 mm lenght
- DTH-1209-365** - 9 mm lenght
- DTH-1210-365** - 10 mm lenght
- DTH-1211-365** - 11 mm lenght
- DTH-1212-365** - 12 mm lenght
- DTH-1213-365** - 13 mm lenght
- DTH-1214-365** - 14 mm lenght
- DTH-1216-365** - 16 mm lenght
- DTH-1218-365** - 18 mm lenght

Thickness: 0,6 mm  
Plates scale: 1/1

Thickness: 0,8 mm  
Plates scale: 1/1


# 1,2 mm



**DTH-0602-41K** Thickness: 0,6mm

**DTH-0802-41K** Thickness: 0,8mm


Avulsion fracture Hook locking plate 1,2mm  
Placa Hook 1,2mm de fratura por avulsão com bloqueio  
Plaque Hook fracture avulsion avec blocage 1,2mm  
Placa Hook 1,2mm para fratura por avulsión con bloqueio



**DTH-0603-41K** Thickness: 0,6mm

**DTH-0803-41K** Thickness: 0,8mm


T Locking Plate 1,2mm - 3 x 8 holes  
Placa 1,2mm em T com bloqueio - 3x8 furos  
Plaque 1,2mm en T avec blocage - 3x8 trous  
Placa de bloqueio en T 1,2mm - 3x8 agujeros



**DTH-0604-41K** Thickness: 0,6mm

**DTH-0804-41K** Thickness: 0,8mm

T Locking Plate 1,2mm - 4 x 8 holes  
Placa 1,2mm em T com bloqueio - 4x8 furos  
Plaque 1,2mm en T avec blocage - 4x8 trous  
Placa de bloqueio en T 1,2mm - 4x8 agujeros



**DTH-0621-41K** Thickness: 0,6mm

**DTH-0821-41K** Thickness: 0,8mm

Y Locking Plate 1,2mm  
Placa 1,2mm em Y com bloqueio  
Plaque 1,2mm en Y avec blocage  
Placa de bloqueio en Y 1,2mm



**DTH-0606-41K** Thickness: 0,6mm

**DTH-0806-41K** Thickness: 0,8mm


Straight Locking Plate 1,2mm - 6 holes  
Placa 1,2mm Reta com Bloqueio - 6 furos  
Plaque 1,2mm droite avec blocage - 6 trous  
Placa de bloqueio recta 1,2 mm - 6 agujeros



**DTH-0612-41K** Thickness: 0,6mm

**DTH-0812-41K** Thickness: 0,8mm


Straight Locking Plate 1,2mm - 12 holes  
Placa 1,2mm Reta com Bloqueio - 12 furos  
Plaque 1,2mm droite avec blocage - 12 trous  
Placa de bloqueio recta 1,2 mm - 12 agujeros



**DTH-0660-41K** Thickness: 0,6mm

**DTH-0860-41K** Thickness: 0,8mm


L Locking Plate 1,2mm -LEFT  
Placa 1,2mm em L com bloqueio - esquerda  
Plaque 1,2mm en L avec blocage - gauche  
Placa de bloqueio en L 1,2mm - izquierda



**DTH-0662-41K** Thickness: 0,6mm

**DTH-0862-41K** Thickness: 0,8mm


L Locking Plate 1,2mm -RIGHT  
Placa 1,2mm em L com bloqueio - direita  
Plaque 1,2mm en L avec blocage - droit  
Placa de bloqueio en L 1,2mm - derecha



**DTH-0670-41K** Thickness: 0,6mm

**DTH-0870-41K** Thickness: 0,8mm

Grid Locking Plate 1,2mm - LEFT  
Placa 1,2mm em rede com bloqueio - esquerda  
Plaque 1,2mm de maille avec blocage - gauche  
Placa de bloqueio malla 1,2mm - izquierda



**DTH-0672-41K** Thickness: 0,6mm

**DTH-0872-41K** Thickness: 0,8mm

Grid Locking Plate 1,2mm - RIGHT  
Placa 1,2mm em rede com bloqueio - direita  
Plaque 1,2mm de maille avec blocage - droit  
Placa de bloqueio malla 1,2mm - derecha

# 1,6 mm

## REGULAR HAND PLATE SYSTEM



compression holes

locking holes

0,6 mm & 0,8 mm thickness plates

 0,6 mm

 0,8 mm



### Self-tapping 1,6 mm screw

Parafuso autorroscante 1,6 mm  
Vis autotaraudante 1,6 mm  
Tornillo autorroscante 1,6 mm



### Self-tapping locking 1,6 mm screw


Parafuso de bloqueio autorroscante 1,6 mm  
Vis de bloqueio autotaraudante 1,6 mm  
Tornillo de bloqueio autorroscante 1,6 mm



### Self-tapping LAG 1,6 mm screw


Parafuso LAG autorroscante 1,6 mm  
Vis LAG autotaraudante 1,6 mm  
Tornillo LAG autorroscante 1,6 mm

atraumatic tip




1,6 mm  
Ø thread

- DTH-1604-265** - 4 mm lenght
- DTH-1605-265** - 5 mm lenght
- DTH-1606-265** - 6 mm lenght
- DTH-1607-265** - 7 mm lenght
- DTH-1608-265** - 8 mm lenght
- DTH-1609-265** - 9 mm lenght
- DTH-1610-265** - 10 mm lenght
- DTH-1611-265** - 11 mm lenght
- DTH-1612-265** - 12 mm lenght
- DTH-1613-265** - 13 mm lenght
- DTH-1614-265** - 14 mm lenght
- DTH-1616-265** - 16 mm lenght
- DTH-1618-265** - 18 mm lenght
- DTH-1620-265** - 20 mm lenght
- DTH-1622-265** - 22 mm lenght
- DTH-1624-265** - 24 mm lenght



1,6 mm  
Ø thread

- DTH-1604-265K** - 4 mm lenght
- DTH-1605-265K** - 5 mm lenght
- DTH-1606-265K** - 6 mm lenght
- DTH-1607-265K** - 7 mm lenght
- DTH-1608-265K** - 8 mm lenght
- DTH-1609-265K** - 9 mm lenght
- DTH-1610-265K** - 10 mm lenght
- DTH-1611-265K** - 11 mm lenght
- DTH-1612-265K** - 12 mm lenght
- DTH-1613-265K** - 13 mm lenght
- DTH-1614-265K** - 14 mm lenght
- DTH-1616-265K** - 16 mm lenght
- DTH-1618-265K** - 18 mm lenght
- DTH-1620-265K** - 20 mm lenght
- DTH-1622-265K** - 22 mm lenght
- DTH-1624-265K** - 24 mm lenght




1,6 mm  
Ø thread

- DTH-1606-365** - 6 mm lenght
- DTH-1607-365** - 7 mm lenght
- DTH-1608-365** - 8 mm lenght
- DTH-1609-365** - 9 mm lenght
- DTH-1610-365** - 10 mm lenght
- DTH-1611-365** - 11 mm lenght
- DTH-1612-365** - 12 mm lenght
- DTH-1613-365** - 13 mm lenght
- DTH-1614-365** - 14 mm lenght
- DTH-1616-365** - 16 mm lenght
- DTH-1618-365** - 18 mm lenght
- DTH-1620-365** - 20 mm lenght
- DTH-1622-365** - 22 mm lenght
- DTH-1624-365** - 24 mm lenght


Thickness: 0,6 mm  
Plates scale: 1/1

Thickness: 0,8 mm  
Plates scale: 1/1

# 1,6 mm



**DTH-0603-51K** Thickness: 0,6mm  
**DTH-0803-51K** Thickness: 0,8mm  
T Locking Plate 1,6mm - 3 x 8 holes  
Placa 1,6mm em T com bloqueio - 3x8 furos  
Plaque 1,6mm en T avec blocage - 3x8 trous  
Placa de bloqueo en T 1,6mm - 3x8 agujeros




**DTH-0605-70K** Thickness: 0,6mm  
**DTH-0805-70K** Thickness: 0,8mm  
Subcondylar Locking Plate 1,6mm  
Placa 1,6mm Subcondilar com bloqueio  
Plaque 1,6mm Subcondyle avec blocage  
Placa de bloqueo subcondilar 1,6 mm




**DTH-0606-51K** Thickness: 0,6mm  
**DTH-0806-51K** Thickness: 0,8mm  
Straight Locking Plate 1,6mm - 6 Holes  
Placa 1,6mm Reta com Bloqueio - 6 furos  
Plaque 1,6mm droite avec blocage - 6 trous  
Placa de bloqueo recta 1,6 mm - 6 agujeros




**DTH-0660-51K** Thickness: 0,6mm  
**DTH-0860-51K** Thickness: 0,8mm  
L Locking Plate 1,6mm - LEFT  
Placa 1,6mm em L com bloqueio - esquerda  
Plaque 1,6mm en L avec blocage - gauche  
Placa de bloqueo en L 1,6mm - izquierda




**DTH-0670-51K** Thickness: 0,6mm  
**DTH-0870-51K** Thickness: 0,8mm  
Grid Locking Plate 1,6mm - LEFT  
Placa 1,6mm em rede com bloqueio - esquerda  
Plaque 1,6mm de maille avec blocage - gauche  
Placa de bloqueo malla 1,6mm - izquierda



**DTH-0621-31K** Thickness: 0,6mm  
**DTH-0821-31K** Thickness: 0,8mm  
Condylar Locking Plate 1,6mm - LEFT  
Placa 1,6mm Condilar com bloqueio - esquerda  
Plaque 1,6mm Condyle avec blocage - gauche  
Placa de bloqueo condilar 1,6 mm - izquierda



**DTH-0604-51K** Thickness: 0,6mm  
**DTH-0804-51K** Thickness: 0,8mm  
T Locking Plate 1,6mm - 4 x 8 holes  
Placa 1,6mm em T com bloqueio - 4x8 furos  
Plaque 1,6mm en T avec blocage - 4x8 trous  
Placa de bloqueo en T 1,6mm - 4x8 agujeros




**DTH-0621-51K** Thickness: 0,6mm  
**DTH-0821-51K** Thickness: 0,8mm  
Y Locking Plate 1,6mm  
Placa 1,6mm em Y com bloqueio  
Plaque 1,6mm en Y avec blocage  
Placa de bloqueo en Y 1,6mm




**DTH-0612-51K** Thickness: 0,6mm  
**DTH-0812-51K** Thickness: 0,8mm  
Straight Locking Plate 1,6mm - 12 Holes  
Placa 1,6mm Reta com Bloqueio - 12 furos  
Plaque 1,6mm droite avec blocage - 12 trous  
Placa de bloqueo recta 1,6 mm - 12 agujeros






**DTH-0662-51K** Thickness: 0,6mm  
**DTH-0862-51K** Thickness: 0,8mm  
L Locking Plate 1,6mm - RIGHT  
Placa 1,6mm em L com bloqueio - direita  
Plaque 1,6mm en L avec blocage - droit  
Placa de bloqueo en L 1,6mm - derecha



**DTH-0672-51K** Thickness: 0,6mm  
**DTH-0872-51K** Thickness: 0,8mm  
Grid Locking Plate 1,6mm - RIGHT  
Placa 1,6mm em rede com bloqueio - direita  
Plaque 1,6mm de maille avec blocage - droit  
Placa de bloqueo malla 1,6mm - derecha



**DTH-0623-31K** Thickness: 0,6mm  
**DTH-0823-31K** Thickness: 0,8mm  
Condylar Locking Plate 1,6mm - RIGHT  
Placa 1,6mm Condilar com bloqueio - direita  
Plaque 1,6mm Condyle avec blocage - droit  
Placa de bloqueo condilar 1,6 mm - derecha




**DTH-0602-51K** Thickness: 0,6mm  
**DTH-0802-51K** Thickness: 0,8mm  
Avulsion fracture Hook locking plate 1,6mm  
Placa Hook 1,6mm de fratura por avulsão com bloqueio  
Plaque Hook fracture avulsion avec blocage 1,6mm  
Placa Hook 1,6mm para fratura por avulsión con bloqueio



# 2,0 mm

## REGULAR HAND PLATE SYSTEM



compression holes

locking holes

1,0 mm thickness plates



### Self-tapping 2,0 mm screw

Parafuso autorroscante 2,0 mm  
Vis autotaraudante 2,0 mm  
Tornillo autorroscante 2,0 mm



### Self-tapping locking 2,0 mm screw


Parafuso de bloqueio autorroscante 2,0 mm  
Vis de bloqueio autotaraudante 2,0 mm  
Tornillo de bloqueio autorroscante 2,0 mm



### Self-tapping LAG 2,0 mm screw


Parafuso LAG autorroscante 2,0 mm  
Vis LAG autotaraudante 2,0 mm  
Tornillo LAG autorroscante 2,0 mm

atraumatic tip




2,0 mm  
Ø thread

<b>DTH-2006-265</b>	- 6 mm lenght
<b>DTH-2007-265</b>	- 7 mm lenght
<b>DTH-2008-265</b>	- 8 mm lenght
<b>DTH-2009-265</b>	- 9 mm lenght
<b>DTH-2010-265</b>	- 10 mm lenght
<b>DTH-2011-265</b>	- 11 mm lenght
<b>DTH-2012-265</b>	- 12 mm lenght
<b>DTH-2013-265</b>	- 13 mm lenght
<b>DTH-2014-265</b>	- 14 mm lenght
<b>DTH-2016-265</b>	- 16 mm lenght
<b>DTH-2018-265</b>	- 18 mm lenght
<b>DTH-2020-265</b>	- 20 mm lenght
<b>DTH-2022-265</b>	- 22 mm lenght
<b>DTH-2024-265</b>	- 24 mm lenght
<b>DTH-2026-265</b>	- 26 mm lenght
<b>DTH-2028-265</b>	- 28 mm lenght
<b>DTH-2030-265</b>	- 30 mm lenght
<b>DTH-2032-265</b>	- 32 mm lenght
<b>DTH-2034-265</b>	- 34 mm lenght
<b>DTH-2036-265</b>	- 36 mm lenght



2,0 mm  
Ø thread

<b>DTH-2006-265K</b>	- 6 mm lenght
<b>DTH-2007-265K</b>	- 7 mm lenght
<b>DTH-2008-265K</b>	- 8 mm lenght
<b>DTH-2009-265K</b>	- 9 mm lenght
<b>DTH-2010-265K</b>	- 10 mm lenght
<b>DTH-2011-265K</b>	- 11 mm lenght
<b>DTH-2012-265K</b>	- 12 mm lenght
<b>DTH-2013-265K</b>	- 13 mm lenght
<b>DTH-2014-265K</b>	- 14 mm lenght
<b>DTH-2016-265K</b>	- 16 mm lenght
<b>DTH-2018-265K</b>	- 18 mm lenght
<b>DTH-2020-265K</b>	- 20 mm lenght
<b>DTH-2022-265K</b>	- 22 mm lenght
<b>DTH-2024-265K</b>	- 24 mm lenght
<b>DTH-2026-265K</b>	- 26 mm lenght
<b>DTH-2028-265K</b>	- 28 mm lenght
<b>DTH-2030-265K</b>	- 30 mm lenght
<b>DTH-2032-265K</b>	- 32 mm lenght
<b>DTH-2034-265K</b>	- 34 mm lenght
<b>DTH-2036-265K</b>	- 36 mm lenght




2,0 mm  
Ø thread

<b>DTH-2006-365</b>	- 6 mm lenght
<b>DTH-2007-365</b>	- 7 mm lenght
<b>DTH-2008-365</b>	- 8 mm lenght
<b>DTH-2009-365</b>	- 9 mm lenght
<b>DTH-2010-365</b>	- 10 mm lenght
<b>DTH-2011-365</b>	- 11 mm lenght
<b>DTH-2012-365</b>	- 12 mm lenght
<b>DTH-2013-365</b>	- 13 mm lenght
<b>DTH-2014-365</b>	- 14 mm lenght
<b>DTH-2016-365</b>	- 16 mm lenght
<b>DTH-2018-365</b>	- 18 mm lenght
<b>DTH-2020-365</b>	- 20 mm lenght
<b>DTH-2022-365</b>	- 22 mm lenght
<b>DTH-2024-365</b>	- 24 mm lenght
<b>DTH-2026-365</b>	- 26 mm lenght
<b>DTH-2028-365</b>	- 28 mm lenght
<b>DTH-2030-365</b>	- 30 mm lenght
<b>DTH-2032-365</b>	- 32 mm lenght
<b>DTH-2034-365</b>	- 34 mm lenght
<b>DTH-2036-365</b>	- 36 mm lenght


Thickness: 1,0 mm  
Plates scale: 1/1

# 2,0mm





**DTH-1021-61K**

Y Locking Plate 2,0mm  
Placa 2,0mm em Y com bloqueio  
Plaque 2,0mm en Y avec blocage  
Placa de bloqueo en Y 2,0mm




**DTH-1003-61K**

T Locking Plate 2,0mm 3 x 8 holes  
Placa 2,0mm em T com bloqueio - 3x8 furos  
Plaque 2,0mm en T avec blocage - 3x8 trous  
Placa de bloqueo en T 2,0mm - 3x8 agujeros




**DTH-1005-70K**

Subcondylar Locking Plate 2,0mm  
Placa 2,0mm Subcondilar com bloqueio  
Plaque 2,0mm Subcondyle avec blocage  
Placa de bloqueo subcondilar 2,0 mm




**DTH-1002-61K**

T Locking Plate 2,0mm - 2 x 8 holes  
Placa 2,0mm em T com bloqueio - 2x8 furos  
Plaque 2,0mm en T avec blocage - 2x8 trous  
Placa de bloqueo en T 2,0mm - 2x8 agujeros




oblong compression holes  
only for cortical screws




**DTH-1006-61K**

Straight Locking Plate 2,0mm - 6 Holes  
Placa 2,0mm Reta com Bloqueio - 6 furos  
Plaque 2,0mm droite avec blocage - 6 trous  
Placa de bloqueo recta 2,0 mm - 6 agujeros




**DTH-1012-61K**

Straight Locking Plate 2,0mm - 12 Holes  
Placa 2,0mm Reta com Bloqueio - 12 furos  
Plaque 2,0mm droite avec blocage - 12 trous  
Placa de bloqueo recta 2,0 mm - 12 agujeros




**DTH-1060-61K**

L Locking Plate 2,0mm - **LEFT**  
Placa 2,0mm em L com bloqueio - **esquerda**  
Plaque 2,0mm en L avec blocage - **gauche**  
Placa de bloqueo en L 2,0mm - **izquierda**




**DTH-1062-61K**

L Locking Plate 2,0mm - **RIGHT**  
Placa 2,0mm em L com bloqueio - **direita**  
Plaque 2,0mm en L avec blocage - **droit**  
Placa de bloqueo en L 2,0mm - **derecha**




**DTH-1021-31K**

Condylar Locking Plate 2,0mm - **LEFT**  
Placa 2,0mm Condilar com bloqueio - **esquerda**  
Plaque 2,0mm Condyle avec blocage - **gauche**  
Placa de bloqueo condilar 2,0 mm - **izquierda**



**DTH-1023-31K**

Condylar Locking Plate 2,0mm - **RIGHT**  
Placa 2,0mm Condilar com bloqueio - **direita**  
Plaque 2,0mm Condyle avec blocage - **droit**  
Placa de bloqueo condilar 2,0 mm - **derecha**



**DTH-1081-61K**


Z Locking Plate 2,0mm  
Placa 2,0mm em Z com bloqueio  
Plaque 2,0mm en Z avec blocage  
Placa de bloqueo en Z 2,0mm






# 2,0mm

Thickness: 1,0 mm  
Plates scale: 1/1




### DTH-1070-61K

Grid Locking Plate 2,0mm - **LEFT**  
Placa 2,0mm em rede com bloqueio - **esquerda**  
Plaque 2,0mm de maille avec blocage - **gauche**  
Placa de bloqueo malla 2,0mm - **izquierda**




### DTH-1072-61K

Grid Locking Plate 2,0mm - **RIGHT**  
Placa 2,0mm em rede com bloqueio - **direita**  
Plaque 2,0mm de maille avec blocage - **droit**  
Placa de bloqueo malla 2,0mm - **derecha**




### DTH-1004-60K

Straight LCP\* 2,0mm - 4 Holes  
Placa 2,0mm reta LCP\* - 4 furos  
Plaque 2,0mm LCP\* droite - 4 trous  
Placa recta LCP\* 2,0mm - 4 agujeros



### DTH-1006-60K


Straight LCP\* 2,0mm - 6 Holes  
Placa 2,0mm reta LCP\* - 6 furos  
Plaque 2,0mm LCP\* droite - 6 trous  
Placa recta LCP\* 2,0mm - 6 agujeros



### DTH-1008-60K


Straight LCP\* 2,0mm - 8 Holes  
Placa 2,0mm reta LCP\* - 8 furos  
Plaque 2,0mm LCP\* droite - 8 trous  
Placa recta LCP\* 2,0mm - 8 agujeros

\*LCP: Local Compression Plate



# 2,4mm

## REGULAR HAND PLATE SYSTEM



compression holes

locking holes

1,3 mm thickness plates



### Self-tapping 2,4 mm screw

Parafuso autorroscante 2,4 mm  
Vis autotaraudante 2,4 mm  
Tornillo autorroscante 2,4 mm



### Self-tapping locking 2,4 mm screw


Parafuso de bloqueio autorroscante 2,4 mm  
Vis de blocage autotaraudante 2,4 mm  
Tornillo de bloqueo autorroscante 2,4 mm



atraumatic tip


### Self-tapping LAG 2,4 mm screw

Parafuso LAG autorroscante 2,4 mm  
Vis LAG autotaraudante 2,4 mm  
Tornillo LAG autorroscante 2,4 mm




2,4 mm  
Ø thread

- DTH-2406-265** - 6 mm lenght
- DTH-2407-265** - 7 mm lenght
- DTH-2408-265** - 8 mm lenght
- DTH-2409-265** - 9 mm lenght
- DTH-2410-265** - 10 mm lenght
- DTH-2411-265** - 11 mm lenght
- DTH-2412-265** - 12 mm lenght
- DTH-2413-265** - 13 mm lenght
- DTH-2414-265** - 14 mm lenght
- DTH-2416-265** - 16 mm lenght
- DTH-2418-265** - 18 mm lenght
- DTH-2420-265** - 20 mm lenght
- DTH-2422-265** - 22 mm lenght
- DTH-2424-265** - 24 mm lenght
- DTH-2426-265** - 26 mm lenght
- DTH-2428-265** - 28 mm lenght
- DTH-2430-265** - 30 mm lenght
- DTH-2432-265** - 32 mm lenght
- DTH-2434-265** - 34 mm lenght
- DTH-2436-265** - 36 mm lenght




2,4 mm  
Ø thread

- DTH-2406-265K** - 6 mm lenght
- DTH-2407-265K** - 7 mm lenght
- DTH-2408-265K** - 8 mm lenght
- DTH-2409-265K** - 9 mm lenght
- DTH-2410-265K** - 10 mm lenght
- DTH-2411-265K** - 11 mm lenght
- DTH-2412-265K** - 12 mm lenght
- DTH-2413-265K** - 13 mm lenght
- DTH-2414-265K** - 14 mm lenght
- DTH-2416-265K** - 16 mm lenght
- DTH-2418-265K** - 18 mm lenght
- DTH-2420-265K** - 20 mm lenght
- DTH-2422-265K** - 22 mm lenght
- DTH-2424-265K** - 24 mm lenght
- DTH-2426-265K** - 26 mm lenght
- DTH-2428-265K** - 28 mm lenght
- DTH-2430-265K** - 30 mm lenght
- DTH-2432-265K** - 32 mm lenght
- DTH-2434-265K** - 34 mm lenght
- DTH-2436-265K** - 36 mm lenght




2,4 mm  
Ø thread

- DTH-2406-365** - 6 mm lenght
- DTH-2407-365** - 7 mm lenght
- DTH-2408-365** - 8 mm lenght
- DTH-2409-365** - 9 mm lenght
- DTH-2410-365** - 10 mm lenght
- DTH-2411-365** - 11 mm lenght
- DTH-2412-365** - 12 mm lenght
- DTH-2413-365** - 13 mm lenght
- DTH-2414-365** - 14 mm lenght
- DTH-2416-365** - 16 mm lenght
- DTH-2418-365** - 18 mm lenght
- DTH-2420-365** - 20 mm lenght
- DTH-2422-365** - 22 mm lenght
- DTH-2424-365** - 24 mm lenght
- DTH-2426-365** - 26 mm lenght
- DTH-2428-365** - 28 mm lenght
- DTH-2430-365** - 30 mm lenght
- DTH-2432-365** - 32 mm lenght
- DTH-2434-365** - 34 mm lenght
- DTH-2436-365** - 36 mm lenght




Thickness: 1,3 mm  
Plates scale: 1/1

# 2,4mm




**DTH-1321-60K**

Y Locking Plate 2,4mm  
Placa 2,4mm em Y com bloqueio  
Plaque 2,4mm en Y avec blocage  
Placa de bloqueo en Y 2,4mm




**DTH-1303-61K**

T Locking Plate 2,4mm - 3 x 8 holes  
Placa 2,4mm em T com bloqueio - 3x8 furos  
Plaque 2,4mm en T avec blocage - 3x8 trous  
Placa de bloqueo en T 2,4mm - 3x8 agujeros




**DTH-1305-70K**

Subcondylar Locking Plate 2,4mm  
Placa 2,4mm Subcondilar com bloqueio  
Plaque 2,4mm Subcondyle avec blocage  
Placa de bloqueo subcondilar 2,4mm




**DTH-1302-61K**

T Locking Plate 2,4mm - 2 x 8 holes  
Placa 2,4mm em T com bloqueio - 2x8 furos  
Plaque 2,4mm en T avec blocage - 2x8 trous  
Placa de bloqueo en T 2,4mm - 2x8 agujeros




**DTH-1306-61K**

Straight Locking Plate 2,4mm - 6 holes  
Placa 2,4mm Reta com bloqueio - 6 furos  
Plaque 2,4mm droite avec blocage - 6 trous  
Placa de bloqueo recta 2,4 mm - 6 agujeros




**DTH-1312-61K**

Straight Locking Plate 2,4mm - 12 holes  
Placa 2,4 mm Reta com bloqueio - 12 furos  
Plaque 2,4mm droite avec blocage - 12 trous  
Placa de bloqueo recta 2,4 mm - 12 agujeros




**DTH-1360-61K**

L Locking Plate 2,4mm - LEFT  
Placa 2,4mm em L com bloqueio - esquerda  
Plaque 2,4mm en L avec blocage - gauche  
Placa de bloqueo en L 2,4mm - izquierda




**DTH-1362-61K**

L Locking Plate 2,4mm - RIGHT  
Placa 2,4mm em L com bloqueio - direita  
Plaque 2,4mm en L avec blocage - droit  
Placa de bloqueo en L 2,4mm - derecha




**DTH-1321-61K**

Condylar Locking Plate 2,4mm - LEFT  
Placa 2,4mm Condilar com bloqueio - esquerda  
Plaque 2,4mm Condyle avec blocage - gauche  
Placa de bloqueo condilar 2,4 mm - izquierda



**DTH-1323-61K**

Condylar Locking Plate 2,4mm - RIGHT  
Placa 2,4mm Condilar com bloqueio - direita  
Plaque 2,4mm Condyle avec blocage - droit  
Placa de bloqueo condilar 2,4 mm - derecha




**DTH-1381-61K**

Z Locking Plate 2,4mm  
Placa 2,4mm em Z com bloqueio  
Plaque 2,4mm en Z avec blocage  
Placa de bloqueo en Z 2,4mm




# 2,4mm

Thickness: 1,3 mm  
Plates scale: 1/1




### DTH-1370-61K

Grid Locking Plate 2,4mm - **LEFT**  
Placa 2,4mm em rede com bloqueio - **esquerda**  
Plaque 2,4mm de maille avec blocage - **gauche**  
Placa de bloqueo malla 2,4mm - **izquierda**




### DTH-1372-61K

Grid Locking Plate 2,4mm - **RIGHT**  
Placa 2,4mm em rede com bloqueio - **direita**  
Plaque 2,4mm de maille avec blocage - **droit**  
Placa de bloqueo malla 2,4mm - **derecha**




### DTH-1304-60K

Straight LCP\* 2,4mm - 4 Holes  
Placa 2,4mm reta LCP\* - 4 furos  
Plaque 2,4mm LCP\* droite - 4 trous  
Placa recta LCP\* 2,4mm - 4 agujeros



### DTH-1306-60K


Straight LCP\* 2,4mm - 6 Holes  
Placa 2,4mm reta LCP\* - 6 furos  
Plaque 2,4mm LCP\* droite - 6 trous  
Placa recta LCP\* 2,4mm - 6 agujeros



### DTH-1308-60K

Straight LCP\* 2,4mm - 8 Holes  
Placa 2,4mm reta LCP\* - 8 furos  
Plaque 2,4mm LCP\* droite - 8 trous  
Placa recta LCP\* 2,4mm - 8 agujeros

\*LCP: Local Compression Plate





**DELPHOS**  
Implants

*Embracing Life™*

**ATTIS**  
MEDICAL